

KD-Validated Anti-ALPL Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1924

Product Overview

Name	KD-Validated Anti-ALPL Mouse Monoclonal Antibody
Accession(Primary)	P05186
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ALPL
Species Reactivity	Human
Calculated MW	Predicted, 57 kDa, observed, 75 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P05186)

Gene Name	ALPL
About Gene Symbol	ALPL (Alkaline Phosphatase, Biom mineralization Associated) is a protein phosphatase involved in cellular metabolism. It is an important research target in bone biology and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

