

KD-Validated Anti-PPP5C Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1926

Product Overview


Name	KD-Validated Anti-PPP5C Mouse Monoclonal Antibody
Accession(Primary)	P53041
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Immunogen	Recombinant protein of human PPP5C
Species Reactivity	Human
Calculated MW	Predicted, 57 kDa, observed, 52 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:2,500

Target information (P53041)

Gene Name	PPP5C
About Gene Symbol	PPP5C (Protein Phosphatase 5 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024, GenScript Biotechnologies LLC



Copyright ©2024, GenScript Biotechnologies LLC