

KD-Validated Anti-G protein subunit beta 2 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1932

Product Overview

Name	KD-Validated Anti-G protein subunit beta 2 Rabbit Monoclonal Antibody
Accession(Primary)	P62879
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GNB2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 37 kDa , observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P62879)

Gene Name	GNB2
About Gene Symbol	GNB2 (Transducin Beta Chain 2) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellu

