

KD-Validated Anti-APMAP Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1939

Product Overview

Name	KD-Validated Anti-APMAP Mouse Monoclonal Antibody
Accession(Primary)	Q9HDC9
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human C20orf3
Species Reactivity	Human
Calculated MW	Predicted, 46 kDa, observed, 46 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000


Target information (Q9HDC9)

Gene Name	APMAP
About Gene Symbol	APMAP (Adipocyte Plasma Membrane Associated Protein) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-APMAP antibody (Cat#ABG1939). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-APMAP antibody (Cat#ABG1939, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-APMAP antibody (Cat#ABG1939). APMAP expression in wild type (WT) and APMAP shRNA knockdown (KD) HT-1080 cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-APMAP antibody (Cat#ABG1939, 1:5,000) and HRP-conjugated goat anti-mouse secondary antibody respectively.