

KD-Validated Anti-Chromogranin A Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1943

Product Overview

Name	KD-Validated Anti-Chromogranin A Mouse Monoclonal Antibody
Accession(Primary)	P10645
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human CHGA
Species Reactivity	Human
Calculated MW	Predicted, 51 kDa, observed, 48 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000Immunocytochemistry (IC): 1:100-1:1,000

Target information (P10645)

Gene Name	CHGA
About Gene Symbol	CHGA (Chromogranin A) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

