

KD-Validated Anti-SERPINB1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1944

Product Overview


Name	KD-Validated Anti-SERPINB1 Mouse Monoclonal Antibody
Accession(Primary)	P30740
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Immunogen	Recombinant protein of human SERPINB1
Species Reactivity	Human
Calculated MW	Predicted, 43 kDa, observed, 43 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P30740)


Gene Name	SERPINB1
About Gene Symbol	SERPINB1 (Serp Family B Member 1) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC

Serpin family B member 1-Alexa Fluor® 647