

KD-Validated Anti-ETFA Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1952

Product Overview


Name	KD-Validated Anti-ETFA Mouse Monoclonal Antibody
Accession(Primary)	P13804
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ETFA
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 35 kDa, observed, 30 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:2,500-1:5,000

Target information (P13804)

Gene Name	ETFA
About Gene Symbol	ETFA (Electron Transfer Flavoprotein Subunit Alpha) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2014 Genium Biotechnologies LLC



Copyright © 2014 Genium Biotechnologies LLC