

KD-Validated Anti-HADHA Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1958

Product Overview

Name	KD-Validated Anti-HADHA Mouse Monoclonal Antibody
Accession(Primary)	P40939
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human HADHA
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 83 kDa, observed, 73 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:2,500 Flow cytometry (FCM): 1:2,000


Target information (P40939)

Gene Name	HADHA
About Gene Symbol	HADHA is a mitochondrial acyltransferase enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic




Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-HADHA antibody (Cat#ABG1958). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HADHA antibody (Cat#ABG1958, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-HADHA antibody (Cat#ABG1958). HADHA expression in wild type (WT) and HADHA shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-HADHA antibody (Cat#ABG1958, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Flow cytometric analysis of GIPC1 expression in HepG2 cells using anti-GIPC1 antibody (Cat#63585, 1:2,000). Green, isotype control; red, GIPC1.