

KD-Validated Anti-G protein subunit alpha q Rabbit

Rabbit monoclonal antibody

ABG1960

Product Overview

Name	KD-Validated Anti-G protein subunit alpha q Rabbit Monoclonal Antibody
Accession(Primary)	P50148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GNAQ
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 42 kDa , observed, 42 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P50148)

Gene Name	GNAQ
About Gene Symbol	GNAQ (G Protein Subunit Alpha Q) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cell

