

KD-Validated Anti-PPP5C mouse monoclonal antibody

Mouse monoclonal antibody




ABG1961

Product Overview

Name	KD-Validated Anti-PPP5C mouse monoclonal antibody
Accession(Primary)	P53041
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human PPP5C
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 57 kDa, observed, 55 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P53041)

Gene Name	PPP5C
About Gene Symbol	PPP5C (Protein Phosphatase 5 Catalytic Subunit) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Protein phosphatase 5 catalytic subunit-Alexa Fluor® 647