

KD-Validated Anti-AK3 mouse monoclonal antibody

Mouse monoclonal antibody


ABG1963

Product Overview


Name	KD-Validated Anti-AK3 mouse monoclonal antibody
Accession(Primary)	Q9UIJ7
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human AK3
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 26 kDa, observed, 26 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (Q9UIJ7)

Gene Name	AK3
About Gene Symbol	AK3 (Adenylate Kinase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC