

KD-Validated Anti-ATIC Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1969

Product Overview

Name	KD-Validated Anti-ATIC Mouse Monoclonal Antibody
Accession(Primary)	P31939
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ATIC
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 65 kDa, observed, 65 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P31939)

Gene Name	ATIC
About Gene Symbol	ATIC is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

