

KD-Validated Anti-AXL Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG1979

Product Overview

Name	KD-Validated Anti-AXL Rabbit Monoclonal Antibody
Accession(Primary)	P30530
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	Recombinant protein of human Axl
Species Reactivity	Human, Rat
Calculated MW	Predicted, 98 kDa, observed, 138 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (P30530)

Gene Name	AXL
About Gene Symbol	AXL (AXL Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

