

KD-Validated Anti-BCR Rabbit Monoclonal Antibody

Rabbit monoclonal antibody




ABG1983


Product Overview

Name	KD-Validated Anti-BCR Rabbit Monoclonal Antibody
Accession(Primary)	P11274
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Bcr
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 143 kDa, observed, 160 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000Immunocytochemistry (IC): 1:100-1:1,000

Target information (P11274)

Gene Name	BCR
About Gene Symbol	BCR (BCR Activator Of RhoGEF And GTPase) is a GTPase-activating protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HepG2 cells with anti-BCR antibody (Cat#ABG1983, 1:1,000). Nuclei were stained blue with DAPI; BCR was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.