

KD-Validated Anti-BID Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1990

Product Overview

Name	KD-Validated Anti-BID Mouse Monoclonal Antibody
Accession(Primary)	P55957
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human BID
Species Reactivity	Human
Calculated MW	Predicted, 22 kDa, observed, 19 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P55957)

Gene Name	BID
About Gene Symbol	BID (BH3 Interacting Domain Death Agonist) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

