

KD-Validated Anti-Transglutaminase 2 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG1994

Product Overview


Name	KD-Validated Anti-Transglutaminase 2 Mouse Monoclonal Antibody
Accession(Primary)	P21980
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human Transglutaminase 2
Species Reactivity	Human, Rat
Calculated MW	Predicted, 77 kDa, observed, 74 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P21980)


Gene Name	TGM2
About Gene Symbol	TGM2 (Transglutaminase 2) is a protease enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC