

KD-Validated Anti-RAB27 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG1995

Product Overview

Name	KD-Validated Anti-RAB27 Mouse Monoclonal Antibody
Accession(Primary)	P51159
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human RAB27
Species Reactivity	Human
Calculated MW	Predicted, 25 kDa, observed, 27 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P51159)

Gene Name	RAB27A
About Gene Symbol	RAB27A (RAB27A, Member RAS Oncogene Family) is a small GTPase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

