

KD-Validated Anti-Lin-28 Homolog B Mouse Monoclonal Antibody

Mouse monoclonal antibody





ABG1999


Product Overview

Name	KD-Validated Anti-Lin-28 Homolog B Mouse Monoclonal Antibody
Accession(Primary)	Q6ZN17
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant protein of human LIN28B
Species Reactivity	Human
Calculated MW	Predicted, 27 kDa, observed, 32 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000 Immunocytochemistry (IC): 1:100-1:1,000

Target information (Q6ZN17)

Gene Name	LIN28B
About Gene Symbol	LIN28B (Lin-28 Homolog B) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.





Immunocytochemical staining of HAP1 cells with anti-Lin-28 homolog B antibody (Cat#ABG1999, 1:1,000). Nuclei were stained blue with DAPI; Lin-28 homolog B was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 µm.