

KD-Validated Anti-GAPDH Rabbit Polyclonal Antibody

Rabbit polyclonal antibody

ABG2001

Product Overview

Name	KD-Validated Anti-GAPDH Rabbit Polyclonal Antibody
Accession(Primary)	P04406
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GAPDH
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 36 kDa, observed, 36 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P04406)

Gene Name	GAPDH
About Gene Symbol	GAPDH (Glyceraldehyde-3-Phosphate Dehydrogenase) is a metabolic dehydrogenase involved in gene transcription and regulation. It is an important research target in aging research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

