

# KD-Validated Anti-Tumor Protein P53 Rabbit Polyclonal

*Rabbit polyclonal antibody*


ABG2002

## Product Overview


|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | KD-Validated Anti-Tumor Protein P53 Rabbit Polyclonal Antibody            |
| <b>Accession(Primary)</b>   | P04637                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Polyclonal                                                                |
| <b>Isotype</b>              | Rabbit IgG                                                                |
| <b>Immunogen</b>            | A synthesized peptide derived from human p53                              |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Predicted, 4 kDa, observed, 52 kDa                                        |
| <b>Storage Buffer</b>       | Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide. |
| <b>Recommended Dilution</b> | Western Blotting (WB): 1:1,000-1:5,000                                    |

## Target information (P04637)

|                          |                                                                                                                                                                                                                                      |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TP53                                                                                                                                                                                                                                 |
| <b>About Gene Symbol</b> | TP53 (Tumor Protein P53) is a tumor suppressor protein involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Copyright ©2024 Genium Biotechnologies LLC



Copyright ©2024 Genium Biotechnologies LLC