

KD-Validated Anti-XIAP Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2004

Product Overview

Name	KD-Validated Anti-XIAP Mouse Monoclonal Antibody
Accession(Primary)	P98170
Host Species	Mouse
Clonality	Monoclonal
Immunogen	Recombinant human XIAP
Species Reactivity	Human
Calculated MW	Predicted, 57 kDa, observed, 49 kDa
Storage Buffer	Supplied in TBS (pH 7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin 300.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P98170)

Gene Name	XIAP
About Gene Symbol	XIAP (X-Linked Inhibitor Of Apoptosis) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway and apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

