

KD-Validated Anti-CNP Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2005

Product Overview

Name	KD-Validated Anti-CNP Mouse Monoclonal Antibody
Accession(Primary)	P09543
Host Species	Mouse
Clonality	Monoclonal
Immunogen	KLH conjugated synthetic peptide derived from human CNPase
Species Reactivity	Human
Calculated MW	Predicted, 48 kDa, observed, 40 kDa
Storage Buffer	Supplied in TBS (pH 7.4) containing 50% glycerol, 1% BSA and 0.02% Proclin 300.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000

Target information (P09543)

Gene Name	CNP
About Gene Symbol	CNP (2',3'-Cyclic Nucleotide 3' Phosphodiesterase) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

