

KD-Validated Anti-CBX8 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2014

Product Overview

Name	KD-Validated Anti-CBX8 Mouse Monoclonal Antibody
Accession(Primary)	Q9HC52
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2b
Immunogen	Recombinant protein of human CBX8
Species Reactivity	Human
Calculated MW	Predicted, 43 kDa, observed, 43 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:200-1:1,000 Flow Cytometry (FCM): 1:1,000

Target information (Q9HC52)

Gene Name	CBX8
About Gene Symbol	CBX8 (Chromobox 8) is an epigenetic and nuclear regulatory protein involved in epigenetic regulation. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

