

KD-Validated Anti-CCDC93 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2015

Product Overview


Name	KD-Validated Anti-CCDC93 Mouse Monoclonal Antibody
Accession(Primary)	Q567U6
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human CCDC93
Species Reactivity	Human
Calculated MW	Predicted, 73 kDa, observed, 73 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:1,000

Target information (Q567U6)


Gene Name	CCDC93
About Gene Symbol	CCDC93 (Coiled-Coil Domain Containing 93) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Bio Technologies LLC



Copyright ©2024 Genuin Bio Technologies LLC



Copyright ©2023 Genuin Bio Technologies LLC

Coiled-coil domain containing 93-Alexa Fluor® 647