

KD-Validated Anti-ATP5C1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2016

Product Overview


Name	KD-Validated Anti-ATP5C1 Mouse Monoclonal Antibody
Accession(Primary)	P36542
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ATP5C1
Species Reactivity	Human
Calculated MW	Predicted, 33 kDa, observed, 33 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:1,000

Target information (P36542)


Gene Name	ATP5F1C
About Gene Symbol	ATP5F1C (ATP Synthase F1 Subunit Gamma) is a biosynthetic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genium Biotechnologies LLC



Copyright ©2024 Genium Biotechnologies LLC



Copyright ©2024 Genium Biotechnologies LLC