

KD-Validated Anti-CHMP5 Mouse Monoclonal Antibody

Mouse monoclonal antibody




ABG2018

Product Overview

Name	KD-Validated Anti-CHMP5 Mouse Monoclonal Antibody
Accession(Primary)	Q9NZZ3
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human CHMP5
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 25 kDa, observed, 30 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:200-1:1,000 Flow Cytometry (FCM): 1:1,000

Target information (Q9NZZ3)

Gene Name	CHMP5
About Gene Symbol	CHMP5 (Charged Multivesicular Body Protein 5) is an epigenetic and nuclear regulatory protein involved in chromatin remodeling and epigenetics and apoptosis. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Charged multivesicular body protein 5-Alexa Fluor® 647