

KD-Validated Anti-CTDSP2 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2020

Product Overview

Name	KD-Validated Anti-CTDSP2 Mouse Monoclonal Antibody
Accession(Primary)	O14595
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human CTDSP2
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 31 kDa, observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (O14595)

Gene Name	CTDSP2
About Gene Symbol	CTDSP2 (CTD Small Phosphatase 2) is a protein phosphatase involved in cellular metabolism. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

