

KD-Validated Anti-EIF4E2 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2027

Product Overview


Name	KD-Validated Anti-EIF4E2 Mouse Monoclonal Antibody
Accession(Primary)	O60573
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human EIF4E2
Species Reactivity	Human, Mouse
Calculated MW	Predicted, 28 kDa, observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:1,000

Target information (O60573)


Gene Name	EIF4E2
About Gene Symbol	EIF4E2 (Eukaryotic Translation Initiation Factor 4E Family Member 2) is a protein kinase and signaling regulator involved in translational regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is wid



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2024 Genuin Biotechnologies LLC



Copyright © 2025 Genuin Biotechnologies LLC

