

KD-Validated Anti-Epsin 2 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2032

Product Overview


Name	KD-Validated Anti-Epsin 2 Mouse Monoclonal Antibody
Accession(Primary)	O95208
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human EPN2
Species Reactivity	Human
Calculated MW	Predicted, 68 kDa, observed, 68 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500 Flow Cytometry (FCM): 1:1,000

Target information (O95208)


Gene Name	EPN2
About Gene Symbol	EPN2 (Epsin 2) is an actin regulatory protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC