

KD-Validated Anti-FAF1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2039

Product Overview


Name	KD-Validated Anti-FAF1 Mouse Monoclonal Antibody
Accession(Primary)	Q9UNN5
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human FAF1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 74 kDa, observed, 74 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500Flow Cytometry (FCM): 1:1,000Immunocytochemistry (IC): 1:100-1:500

Target information (Q9UNN5)


Gene Name	FAF1
About Gene Symbol	FAF1 (Fas Associated Factor 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC




Copyright ©2025 Genuin Biotechnologies LLC



Copyright ©2025 Genuin Biotechnologies LLC

Fas associated factor 1-Alexa Fluor® 647



Immunocytochemical staining of HepG2 cells with anti-Fas associated factor 1 antibody (Cat#ABG2039, 1:500). Nuclei were stained blue with DAPI; Fas associated factor 1 was stained magenta with Alexa Fluor® 647. Images were taken using Leica stellaris 5. Protein abundance based on laser Intensity and smart gain: Medium. Scale bar, 20 μ m.