

KD-Validated Anti-ECT2 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2041

Product Overview

Name	KD-Validated Anti-ECT2 Mouse Monoclonal Antibody
Accession(Primary)	Q9H8V3
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human ECT2
Species Reactivity	Human
Calculated MW	Predicted, 104 kDa, observed, 104 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:1,000

Target information (Q9H8V3)

Gene Name	ECT2
About Gene Symbol	ECT2 (Epithelial Cell Transforming 2) is a formin actin nucleator involved in cytoskeletal organization. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

