

KD-Validated Anti-EXOSC1 Mouse Monoclonal Antibody

Mouse monoclonal antibody


ABG2043

Product Overview


Name	KD-Validated Anti-EXOSC1 Mouse Monoclonal Antibody
Accession(Primary)	Q9Y3B2
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human EXOSC1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 21 kDa, observed, 21 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q9Y3B2)

Gene Name	EXOSC1
About Gene Symbol	EXOSC1 (Exosome Component 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC