

KD-Validated Anti-LAMP1/CD107a Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2050

Product Overview


Name	KD-Validated Anti-LAMP1 /CD107a Rabbit Monoclonal Antibody
Accession(Primary)	P11279
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human LAMP1
Species Reactivity	Human
Calculated MW	Predicted, 45 kDa , observed, 120 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Immunocytochemistry (IC): 1:1,000

Target information (P11279)


Gene Name	LAMP1
About Gene Symbol	LAMP1 (Lysosomal Associated Membrane Protein 1) is a cell surface marker involved in adaptive immunity and vesicular trafficking. It is an important research target in immunology, flow cytometry, and cancer immunotherapy and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC



Copyright ©2024 Genuin Biotechnologies LLC