

KD-Validated Anti-PI3 Kinase p110 beta Rabbit Monoclonal Antibody

Rabbit monoclonal antibody

ABG2051

Product Overview

Name	KD-Validated Anti-PI3 Kinase p110 beta Rabbit Monoclonal Antibody
Accession(Primary)	P42338
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human PI3 Kinase p110 beta
Species Reactivity	Human
Calculated MW	Predicted , 123 kDa , observed , 110 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000 Flow Cytometry (FCM): 1:2,000

Target information (P42338)

Gene Name	PIK3CB
About Gene Symbol	PIK3CB (Phosphatidylinositol-4,5-Bisphosphate 3-Kinase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely s

