

KD-Validated Anti-Glutaredoxin 3 Mouse Monoclonal Antibody

Mouse monoclonal antibody

ABG2069

Product Overview

Name	KD-Validated Anti-Glutaredoxin 3 Mouse Monoclonal Antibody
Accession(Primary)	O76003
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human GLRX3
Species Reactivity	Human
Calculated MW	Predicted, 37 kDa, observed, 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (O76003)

Gene Name	GLRX3
About Gene Symbol	GLRX3 (Glutaredoxin 3) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

