

KD-Validated Anti-HNRNPUL1 Rabbit Monoclonal Antibody

Rabbit monoclonal antibody


ABG2101

Product Overview


Name	KD-Validated Anti-HNRNPUL1 Rabbit Monoclonal Antibody
Accession(Primary)	Q9BUJ2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human E1B-AP5
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 96 kDa, observed, 90-125 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:1,000-1:5,000Flow Cytometry (FCM): 1:2,000

Target information (Q9BUJ2)


Gene Name	HNRNPUL1
About Gene Symbol	HNRNPUL1 (Heterogeneous Nuclear Ribonucleoprotein U Like 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-HNRNPUL1 antibody (Cat#ABG2101). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HNRNPUL1 antibody (Cat#ABG2101, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-HNRNPUL1 antibody (Cat#ABG2101). HNRNPUL1 expression in wild-type (WT) and HNRNPUL1 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-HNRNPUL1 antibody (Cat#ABG2101, 1:5,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Flow cytometric analysis of HNRNPUL1 expression in HAP-1 cells using anti-HNRNPUL1 antibody (Cat#ABG2101, 1:2,000). Green, isotype control; red, HNRNPUL1.