

KD-Validated Anti-Arginase 2 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG2115

Product Overview


Name	KD-Validated Anti-Arginase 2 Rabbit Polyclonal Antibody
Accession(Primary)	P78540
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Arginase 2
Species Reactivity	Human
Calculated MW	Predicted, 39 kDa, observed, 39 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P78540)

Gene Name	ARG2
About Gene Symbol	ARG2 (Arginase 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-arginase 2 antibody (Cat#ABG2115). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-arginase 2 antibody (Cat#ABG2115, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-arginase 2 antibody (Cat #ABG2115). Arginase 2 expression in wild-type (WT) and arginase 2 (ARG2) shRNA knockdown (KD) HT-1080 cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-arginase 2 antibody (Cat #ABG2115, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody (Cat #201, 1:20,000) respectively.