

# KD-Validated Anti-C1QBP Rabbit Polyclonal Antibody

*Rabbit polyclonal antibody*

ABG2116

## Product Overview

|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | KD-Validated Anti-C1QBP Rabbit Polyclonal Antibody                        |
| <b>Accession(Primary)</b>   | Q07021                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Polyclonal                                                                |
| <b>Isotype</b>              | Rabbit IgG                                                                |
| <b>Immunogen</b>            | A synthesized peptide derived from human BCL2                             |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Predicted, 31 kDa, observed, 35 kDa                                       |
| <b>Storage Buffer</b>       | Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide. |
| <b>Recommended Dilution</b> | Western Blotting (WB): 1:500-1:2,500                                      |

## Target information (Q07021)

|                          |                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | C1QBP                                                                                                                                                                                                                                                                                                                   |
| <b>About Gene Symbol</b> | C1QBP (Complement C1q Binding Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |

