

KD-Validated Anti-Cathepsin D LC Rabbit Polyclonal

Rabbit polyclonal antibody

ABG2120

Product Overview

Name	KD-Validated Anti-Cathepsin D LC Rabbit Polyclonal Antibody
Accession(Primary)	P07339
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cathepsin D LC
Species Reactivity	Human
Calculated MW	Predicted, 45 kDa, observed, 14 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P07339)

Gene Name	CTSD
About Gene Symbol	CTSD (Cathepsin D) is a lysosomal cathepsin protease involved in vesicular trafficking. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

