

KD-Validated Anti-GOT1 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG2121

Product Overview

Name	KD-Validated Anti-GOT1 Rabbit Polyclonal Antibody
Accession(Primary)	P17174
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human GOT1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 46 kDa, observed, 41 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500


Target information (P17174)

Gene Name	GOT1
About Gene Symbol	GOT1 (Glutamic-Oxaloacetic Transaminase 1) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and ide



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-GOT1 antibody (Cat #ABG2121). GOT1 expression in wild-type (WT) and GOT1 shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-GOT1 antibody (Cat #ABG2121, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright © 2024 Genuin Biotechnologies LLC

Western blotting analysis using anti-GOT1 antibody (Cat#ABG2121). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-GOT1 antibody (Cat#ABG2121, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.