

KD-Validated Anti-Galectin 1 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody


ABG2126

Product Overview


Name	KD-Validated Anti-Galectin 1 Rabbit Polyclonal Antibody
Accession(Primary)	P09382
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Galectin 1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 15 kDa, observed, 15 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P09382)

Gene Name	LGALS1
About Gene Symbol	LGALS1 (Galectin 1) is an extracellular matrix laminin protein involved in MAPK/ERK signaling. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-galectin 1 antibody (Cat#ABG2126). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-galectin 1 antibody (Cat#ABG2126, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Western blotting analysis using anti-galectin 1 antibody (Cat #ABG2126). Galectin 1 expression in wild-type (WT) and galectin 1 (LGALS1) shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-galectin 1 antibody (Cat #ABG2126, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.