

KD-Validated Anti-SRC Rabbit Polyclonal Antibody

Rabbit polyclonal antibody

ABG2146

Product Overview

Name	KD-Validated Anti-SRC Rabbit Polyclonal Antibody
Accession(Primary)	P12931
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human c-SRC
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 60 kDa, observed, 58 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P12931)

Gene Name	SRC
About Gene Symbol	SRC (SRC Proto-Oncogene, Non-Receptor Tyrosine Kinase) is a receptor tyrosine kinase involved in receptor-mediated signaling. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

