

KD-Validated Anti-MAPK14 Rabbit Polyclonal Antibody

Rabbit polyclonal antibody

ABG2153

Product Overview

Name	KD-Validated Anti-MAPK14 Rabbit Polyclonal Antibody
Accession(Primary)	Q16539
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human p38
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 41 kDa, observed, 38 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q16539)

Gene Name	MAPK14
About Gene Symbol	MAPK14 (Mitogen-Activated Protein Kinase 14) is a mitogen-activated protein kinase involved in MAPK/ERK signaling and JAK-STAT/cytokine signaling. It is an important research target in inflammatory disease and is widely studied to elucidate cellular mecha

