

KD-Validated Anti-Cleaved-Caspase 8 Rabbit Polyclonal

Rabbit polyclonal antibody


ABG2156

Product Overview

Name	KD-Validated Anti-Cleaved-Caspase 8 Rabbit Polyclonal Antibody
Accession(Primary)	Q14790
Host Species	Rabbit
Clonality	Polyclonal
Isotype	Rabbit IgG
Immunogen	A synthesized peptide derived from human Cleaved-Caspase 8
Species Reactivity	Human
Calculated MW	Predicted, 55 kDa, observed, 28 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500


Target information (Q14790)

Gene Name	CASP8
About Gene Symbol	CASP8 (Caspase 8) is a cysteine protease (caspase) involved in apoptosis. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-caspase 8 antibody (Cat #ABG2156). Total cell lysates (30 μ g) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-caspase 8 antibody (Cat #ABG2156, 1:2,500) and HRP-conjugated goat anti-rabbit secondary antibody respectively.



Copyright ©2025 Genuin Biotechnologies LLC

Western blotting analysis using anti-caspase 8 antibody (Cat #ABG2156). Caspase 8 expression in wild-type (WT) and caspase 8 (CASP8) shRNA knockdown (KD) HeLa cells with 20 μ g of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-caspase 8 antibody (Cat #ABG2156, 1:1,000) and HRP-conjugated goat anti-rabbit secondary antibody respectively.