

KD-Validated Anti-GSTT2 Mouse Monoclonal Antibody

Mouse monoclonal Antibody

ABG2159

Product Overview

Name	KD-Validated Anti-GSTT2 Mouse Monoclonal Antibody
Accession(Primary)	P0CG29
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human GSTT2
Species Reactivity	Human
Calculated MW	Predicted, 28 kDa, observed, 24 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P0CG29)

Gene Name	GSTT2
About Gene Symbol	GSTT2 (Glutathione S-Transferase Theta 2 (Gene/Pseudogene)) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.

