

KD-Validated Anti-MMAB Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2160

Product Overview


Name	KD-Validated Anti-MMAB Mouse Monoclonal Antibody
Accession(Primary)	Q96EY8
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human MMAB
Species Reactivity	Human
Calculated MW	Predicted, 27 kDa, observed, 23 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q96EY8)

Gene Name	MMAB
About Gene Symbol	MMAB (Metabolism Of Cobalamin Associated B) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and id



Western blotting analysis using anti-MMAB antibody (Cat#ABG2160). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-MMAB antibody (Cat#ABG2160, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Western blotting analysis using anti-MMAB antibody (Cat#ABG2160). MMAB expression in wild-type (WT) and MMAB knockdown (KD) HSHC cells with 20 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-MMAB antibody (Cat#ABG2160, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.