

KD-Validated Anti-PHGDH Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2171

Product Overview


Name	KD-Validated Anti-PHGDH Mouse Monoclonal Antibody
Accession(Primary)	O43175
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human PHGDH
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 57 kDa, observed, 57 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (O43175)

Gene Name	PHGDH
About Gene Symbol	PHGDH (Phosphoglycerate Dehydrogenase) is a metabolic dehydrogenase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2025 Genin Biotechnologies LLC



Copyright © 2025 Genin Biotechnologies LLC