

KD-Validated Anti-CEP55 Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2172

Product Overview


Name	KD-Validated Anti-CEP55 Mouse Monoclonal Antibody
Accession(Primary)	Q53EZ4
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human CEP55
Species Reactivity	Human
Calculated MW	Predicted, 54 kDa, observed, 54 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q53EZ4)

Gene Name	CEP55
About Gene Symbol	CEP55 (Centrosomal Protein 55) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2025 Gemini Biotechnologies LLC



Copyright © 2025 Gemini Biotechnologies LLC