

KD-Validated Anti-PDCL3 Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2173

Product Overview


Name	KD-Validated Anti-PDCL3 Mouse Monoclonal Antibody
Accession(Primary)	Q9H2J4
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG2a
Immunogen	Recombinant protein of human PDCL3
Species Reactivity	Human
Calculated MW	Predicted, 28 kDa, observed, 35 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (Q9H2J4)

Gene Name	PDCL3
About Gene Symbol	PDCL3 (Phosducin Like 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright © 2025 Genuin Biotechnologies LLC



Copyright © 2025 Genuin Biotechnologies LLC