

KD-Validated Anti-HNRNPH1 Mouse Monoclonal Antibody

Mouse monoclonal Antibody


ABG2174

Product Overview


Name	KD-Validated Anti-HNRNPH1 Mouse Monoclonal Antibody
Accession(Primary)	P31943
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1
Immunogen	Recombinant protein of human HNRNPH1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Predicted, 49 kDa, observed, 49 kDa
Storage Buffer	Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.
Recommended Dilution	Western Blotting (WB): 1:500-1:2,500

Target information (P31943)

Gene Name	HNRNPH1
About Gene Symbol	HNRNPH1 (Heterogeneous Nuclear Ribonucleoprotein H1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blotting analysis using anti-HNRNPH1 antibody (Cat#ABG2174). Total cell lysates (30 µg) from various cell lines were loaded and separated by SDS-PAGE. The blot was incubated with anti-HNRNPH1 antibody (Cat#ABG2174, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.



Western blotting analysis using anti-HNRNPH1 antibody (Cat#ABG2174). HNRNPH1 expression in wild-type (WT) and HNRNPH1 shRNA knockdown (KD) HeLa cells with 30 µg of total cell lysates. Hsp90 α serves as a loading control. The blot was incubated with anti-HNRNPH1 antibody (Cat#ABG2174, 1:2,500) and HRP-conjugated goat anti-mouse secondary antibody respectively.